

RadMon for Aces and SciCos

- How it works
- How to respond to automated aborts and alarms
- How to manually abort the TeVatron
- Summary – your responsibilities

The Hardware

- On either side of beampipe @ $z = \pm 4.3$ m sits a Beam Loss Monitor (four total)
- BLM signal (prop. to *dose rate*) amplified and digitized in CAMAC modules in control room
- Signal compared with preset thresholds; if too high, drop 5V enable to TeVatron
- Signal also integrated to keep a record of accumulated dose

Monitoring – ACNET E20

E20

SVX Rad Scaler Readout

◆Pgm_Tools◆

*Global Reset

*Plot FIFO
*Select Display Options
*Display Logged Data
Fifos Recording

	Rate (R/s)	Sum (Rads)
W Inner BLM	0	.0378
W Outer BLM	0	0
E Inner BLM	0	4.192
E Outer BLM	0	3.189

Messages

Welcome to the SVX Loss Monitor Page

Monitoring – Fast Time Plots

- Keep a Fast Time Plot of the **integrated dose** going at all times on the ACNET terminal
- Just click on “SVX” on the E-Z Writer page (E11)
- The four integrated BLM doses are E:SVRAD0 – E:SVRAD3

Alarms and Aborts

- **Automated** aborts triggered by high *dose rates* – 12 rad/s
- **Manual** aborts triggered by high *integrated doses* – thresholds depend on what MCR is currently doing (shots, stores, or studies)
- If more than 18 rad collected in past minute, an **automated integrated dose alarm** will sound – start paying close attention to integrated doses (E20 or FTP) in case you have to manually abort

Manual Aborts

- Three successively higher thresholds (~krad, posted by ACNET terminal) must be crossed before you trigger a manual abort
 1. **CDF Manual Alarm** – SciCo alerts MCR
 2. **MCR Manual Abort** – MCR aborts beam
 3. **CDF Manual Abort** – SciCo aborts beam

...ultimately, safety of SVX is in OUR hands!

How to Manually Abort

- In lieu of crash button, manual aborts effected by *lowering automated abort threshold to below pedestal*
 1. Go to ACNET E48, click on “335”
 2. E:SVBLA1 – change yellow number to 10, click “CAUTION!” to confirm
 3. E:SVBTA1 – change yellow number to 10, click “CAUTION!” to confirm

Recovering from Aborts

- All spelled out in CDF procedure 106 – read it!
 1. Silence sono-alarm
 2. Page RDCO, CDF Ops Manager
 3. Once Ops Manager gives the OK, reset the abort hardware (see web pages)
 4. Make an e-log entry, of course

Summary

- Page RDCO if any kind of alarm, abort, or something seems awry
- If any E:SVRAD* activity during your shift, put an ACNET summary plot in e-log
- Constantly monitor E:SVRAD* whenever beam in TeVatron; during high losses, be checking against the manual alarm/abort thresholds
- SciCo should alert MCR when high losses seen
- Further reading: CDF PROC 106 and the RadMon web pages (linked from Monitoring Global page)