

The New Trigger Inhibit System

Jonathan Lewis
Ace Training
March 26 2004




Why be Inhibited?

- Inhibit system insures good data
 - Stop taking triggers quickly and automatically when a problem is detected
 - Example: HV trip
 - Start as soon as detector is in Good state




Out with the old

- Three crates of NIM and Camac logic
 - Familiar to old guys like me
 - Not so easy to maintain



In with the new



CDF Return XPT card

- New transition module
- 160 inputs from patch panel in 1RR12D
- Flexible!
- Maintained by other people!!!!

Overview

- 160 possible inputs
 - Each can be masked in case of problems that can't be resolved quickly
 - Masking means data are suspect
- Two classes of inputs:
 - Signals from iFIX and slow-control PC
 - Includes trips, but takes seconds(?) to reach inhibit
 - One signal per sub-system
 - Inhibit if not OK for data-taking
 - Fast trip signals from power supplies
 - 1 signal per supply
 - msec time scale
 - minimize chance of corrupt data
 - Most common source of inhibit

Operational Considerations


- HV trips are most common source of inhibit
 - Fast trip comes from power supply
 - PC then detects
 - Supply reset, but PC still holds inhibit during ramp
- Different time scales
 - Caen systems and COT have fast response
 - PSM (Low Voltage) polls many supplies
 - Delayed response
 - Still no software inhibit for CLC, CPR, CCR
 - Need to PAUSE run on trips for these systems
- PC signals are new
 - Used by a few systems.
- Inhibits are for PHYSICS table runs with beam or COSMICS to tape
 - Otherwise use the global disable in Run Control

Control Software

- GUI display above ace console
- One icon for each inhibit input
- Use GUI to control mask
 - Lock out changes when run is active in partition that has booked inhibit VME crate
- Don't use old mask function in iFix!!!
 - This page has been deleted



Inhibit Display



Operation

Partition

- InhibitDisplay program checks if the partition owns b0inh00
- If so, masks inhibits for inputs not associated with components in the partition
 - These masks can't be changed in GUI
- Shift crew can now mask inhibits if we need to run with a
 BAD system

Run

Changes to mask are disabled

End

Changes to mask are allowed

Summary of iFix Signals

Channel	Tag	Computer
1	SVX_HV	svxiicon
2	ISL_HV	svxiicon
3	COT_HV	cot2
4	COT_LV	cot2
5	TOF_HV	tof1
6	TOF_LV	tof1
7	CENTRAL_HV	pisabox
8	CES_HV	muon3
9	CES_LV	voltman
10	XENON_OFF	pisabox
11	PLUG_HV	cdfephv
12	PES_LV	voltman
13	MUON_HV	muon3
14	MUON_LV	voltman
15	CSP_CSX	cdfccu
16	IMU_HV	muon3
17	IMU_LV	voltman
18	CLC_HV	clc
19	RACKS	cdf_s3
20	VME_POWER	voltman
21	CSP_CCU	cdfccu
22	IMU_CCU	cdfccu
23	L00_HV	svxiicon

Response to Problems

- Inhibit is a GIGO system
 - Detector experts responsible for inputs
 - If you get an inhibit, and you don't understand why: Call the detector experts!
 - If you have a problem and you didn't get an inhibit, but you think you should: Call the detector experts!
 - If they think they are sending the correct info but the inhibit is behaving badly, call me.
- If an inhibit is masked, the associated system will be considered BAD for that run.
 - Masking an inhibit is a stop-gap measure while things are repaired.

Inject Humor

Run 1 CDF Anthem

Melody Notre Dame fight song

Fight, fight for old CDF.
What is the matter, have you gone deaf?
Hear the DecTalk calling thee:
"Event timeout on machine three."

Whether the gas-gain be green or red, We will take data till we're dead. While our pedestals go shift of to infinity.